

Health Care & Mobility



Powerchair Range

Portable, tough & compact...Your mobility companion



At Cavendish Health Care, we believe that our products should be of the highest quality and provide our customers with something a little bit different. That's why we've made sure that our powerchairs are lightweight, compact and easily transportable.



Our range of SupaChairs have been specifically designed to help you move around in confined spaces - whether that be in your house, flat or bungalow. The sheer variety of the SupaChair models mean that you can choose one that serves you well both inside and outside of your home.



Don't let the size of these powerchairs fool you. As well as having an acute turning circle and excellent manoeuvrability, each one has been manufactured to provide the user with ample legroom and superior ride quality.



To ensure you feel safe, all SupaChairs are equipped with anti-tip wheels and a highly responsive joystick which will help you maintain control when in use.



The SupaChair Mini has been designed specifically for indoor use. The outstanding manoeuvrability gives you the freedom to move around your home with confidence.

This powerchair provides you with the ability to access those rooms that are limited for space. Comfortably navigate your way around furniture, enter and exit bathrooms, bedrooms and turn 360 degrees within a narrow hallway.







Key Features

- Lightweight fully transportable compact design
- Disassemble in seconds
- Folding foot board
- Anti-tip wheels
- Turning radius of 590mm
- Max user weight 23 stone (150kgs)
- Breaks down into 2 parts
- 2 x 150W motors
- Up to 10 miles travel distance
- Maximum speed 3.2mph
- Slope grade ability up to 6 degrees
- Gear stick control with 360° turning capability
- Battery is easily removed for charging
- 2 battery options available: (SLA) 12v 14ah - 9kg (20lbs) (Lithium) 24v 15ah - 4.0kg (9lbs)

MINI SPECIFICATIONS

Drive wheel size	200mm (8 inches)
Castor wheel size	100mm (4 inches)
Overall Length	850mm (33.5 inches
Overall Width	550mm (21.7 inches
Overall Height	800mm (31.5 inches
Standard Seat Width	400mm (15.8 inches
Standard Seat Depth	380mm (15 inches)
Back Rest Height	320mm (12.6 inches
Controller	Dynamic LiNX LE
Weight Capacity	150kg (330lbs)
Total weight with batteries (SLA)	39.1kg (86lbs)
Weight of motor mount complete	20.5kg (45lbs)
Heaviest Piece	20.5kg (45lbs)

Breaks down in to 2 pieces Weight of seat, including armrest 9.4kg (20.7lbs) Maximum speed SLA battery 3.2mph Turning radius (reduced footboard) 590mm (23.2 inches) Ground clearance 42mm (1.6 inches) Battery (SLA / Lithium) 24V 15 ah Weight of battery (SLA / Lithium) 9.2kg / 4kg (20lbs / 9lbs) Slope grade ability 4 degrees (6 degrees with back anti-tip) Travel distance SLA / Lithium Up to 10km / 16.5km (6.2 miles / 10.2 miles) Battery Charger SLA 2A off board Battery charger Lithium 4A off board Armrest height 200mm [7.9 inches] Dynamic breaking and electromagnetic brake Brake

PAGE 2. PAGE 3.



Safari

Lightweight and compact, it's the perfect powerchair that you'll be able to take almost anywhere.

The SupaChair Safari has a lightweight and compact design which can be assembled or disassembled in a matter of seconds. This powerchair is ideal to be used within your home and for those essential shopping trips, this is where you really feel the benefit from the sprung suspension.







- Complete unit will fit into a car boot with ease
- Assembles/disassemble in seconds
- Folding foot board
- Anti-tip wheels
- Turning radius of 700mm
- Max user weight 23 stone (150kgs)
- Breaks down into 2 parts
- 2 x 150W motors
- Up to 10 miles travel distance
- Variable speed control
- Maximum speed 4mph

- Slope grade ability up to 9 degrees
- Gear stick control with 360° turning capability
- Battery is easily removed for charging
- 2 battery options available: (SLA) 12v 14ah - 9kg (20lbs) (Lithium) 24v 15ah - 4.0kg (9lbs)



The Safari Sport has been designed more for outdoor use. Featuring larger wheels & motor, enhancing its performance.

The Safari Sport is nimble enough to navigate small spaces in your home and with the extra power that the Safari Sport offers, it also gives you the versatility for outdoor use for shopping days and with larger wheels provides extra traction when using on rougher terrain for outdoor excursions.





Key Features

- Lightweight, fully transportable, compact design
- Folding foot board
- Anti-tip wheels
- Turning radius of 730mm
- Max user weight 19 stone
- Up to 10 miles travel distance
- Twin seat columns
- Variable speed control
- Maximum speed 5.6mph
- Breaks down into 2 parts
- Complete unit will fit in most car boots
- Slope grade ability up to 9 degrees
- Gear stick control with 360° turning capability
- Comfortable seat with spring suspension
- Battery is easily removed for charging
- 2 battery options available: (SLA) 12v 14ah - 9kg (20lbs) (Lithium) 24v 15ah - 4.0kg (9lbs)

SPECIFICATIONS

Drive wheel size	200mm (8 inches)
Castor wheel size	178mm (7 inches)
Overall Length	870mm (34.3 inches)
Overall Width	600mm (23.6 inches)
Overall Height	840mm (33 inches)
Standard Seat Width	410mm (16.2 inches)
Standard Seat Depth	400mm (15.8 inches)
Back Rest Height	370mm (14.6 inches)
Controller	Dynamic LiNX LE
Weight Capacity	150kg (330lbs)
Total weight with batteries (SLA)	37.7kg (83lbs)
Weight of motor mount complete	17.9kg (39.4lbs)
Heaviest Piece	17.9kg (39.4lbs)

Breaks down in to 2 pieces Weight of seat, including armrest 10.6kg (23.4lbs) Maximum speed SLA battery 4mph Turning radius (reduced footboard) 700mm [27.6 inches] Ground clearance 70mm (2.76inches) Battery (SLA / Lithium) 24V 15 ah Weight of battery (SLA / Lithium) 9.2kg / 4kg (20lbs / 9lbs) Slope grade ability 9 degrees Up to 10km / 16.5km (6.2 miles / 10.2 miles) Travel distance SLA / Lithium Battery Charger SLA 2A off board Battery charger Lithium 4A off board Armrest height 200mm (7.9 inches) Brake Dynamic breaking and electromagnetic brake

PORT SPECIFICATIONS

(
260mm (10 inches)
200mm (8 inches)
910mm (35.8 inches)
610mm (24 inches)
865mm (34 inches)
410mm (16.2 inches)
400mm (15.8 inches)
370mm (14.6 inches)
Dynamic LiNX LE
125kg (276lbs)
40kg (88lbs)
20.2kg (44lbs)
20.2kg (44lbs)

Breaks down in to	2 pieces
Weight of seat, including armrest	10.6kg (23.4lbs)
Maximum speed SLA battery	5.6mph
Turning radius (reduced footboard)	730mm (28.7 inch
Ground clearance	105 (4.1 inches)
Battery (SLA / Lithium)	24V 15ah
Weight of battery (SLA / Lithium)	9.2kg / 4kg (20lbs
Slope grade ability	8 degrees
Travel distance SLA / Lithium	Up to10km / 16.5k
Battery Charger SLA	2A off board
Battery charger Lithium	4A off board
Armrest height	200mm (7.9 inche
Brake	Dynamic breaking

hesl / 9lbs) ikm (6.2 miles / 10.2 miles) es) ng and electromagnetic brake

PAGE 4. PAGE 5.



The SupaChair COMBI is a manual wheelchair and powerchair in one expertly designed package, which folds away easily for quick storage.

The comfort level and range of adjustments of the Combi exceeds that of most electric wheelchairs and even experienced users are impressed with the performance and all-round capabilities. Available in two versions, the standard and the sport.



Key Features

- Adjustable armrests
- Adjustable slung backrest
- 3 back angle setting
- Folding foot rest
- Anti-tip wheels
- Max user weight 23 stone (150kgs)
- Breaks down into 6 parts

- Swing away joystick
- Up to 10 miles travel distance
- Variable speed control
- Maximum speed 4mph
- Slope grade ability up to 9 degrees
- Gear stick control with 360° turning capability
- Rear mount luggage (optional)
- Battery is easily removed for charging
- 2 battery options available: (SLA) 12v 14ah - 9kg (20lbs) (Lithium) 24v 15ah - 4.0kg (9lbs)

The SupaChair COMBI is so versatile that it can be transformed within seconds into a manual version.

This manual wheelchair is complete with anti-tip wheels at the rear and a high backrest for added comfort. The quick release axles help to quickly fold the wheelchair for storage, perfect for use around the home or longer journeys.





Quick release adjustment





Adjustable for left & right hand use



Joystick for attendant (Optional)







3 angle settings for backrest



Rear mounted luggage frame (Optional)



Detachable foot board (Optional)

COMBI SPECIFICATIONS

Drive wheel size	200mm (8 inc
Castor wheel size	178mm (7 ind
Overall Length	1180mm (46.
Overall Width	620mm (25.2
Overall Height	1010mm (39.
Standard Seat Width	490mm (19.3
Standard Seat Depth	450mm (17.7
Back Rest Height	470mm (18.5
Weight Capacity	150kg (330lbs
Total weight with batteries (SLA)	37.3kg (82lbs
Weight of motor mount complete	15kg (33lbs)
Heaviest Piece	15kg (33lbs)

nm (8 inches) nm (7 inches) mm (46.5 inches) nm (25.2 inches) mm (39.8 inches) nm [19.3 inches] nm (17.7 inches) nm (18.5 inches) (g (330lbs) kg (82lbs)

Breaks down into: Maximum speed SLA battery Turning radius Ground clearance Battery (SLA / Lithium) Weight of battery (SLA / Lithium) Slope grade ability Travel distance SLA / Lithium Battery Charger SLA Battery charger Lithium Armrest height Brake

6 pieces 6.5kph (4mph) 11mm (43.3 inches) 83mm (2.76inches) 24V 15 ah 9.2kg / 4kg (20lbs / 9lbs) 9 degrees Up to 10km / 16.5km (6.2 miles / 10.2 miles) 2A off board 4A off board 160-215mm (6.3 - 8.5 inches) Dynamic breaking and electromagnetic brake

BREAKS DOWN IN 6 EASY TO HANDLE PARTS



Backrest

Battery



Chair Frame



Footrest



Motorised **Drive Mount**

PAGE 6. PAGE 7.



www.supachair.co.uk

Cavendish Healthcare Group

showroom 11 The Parkwood Centre, Aston Road, Waterlooville, Hampshire P07 7HT

call 02392 255 600 email sales@cavendish.healthcare







